



+1 (512) 715-9995
+1 (800) 965-9229 Toll Free
+1 (512) 715-9954 Fax
+1 (800) 806-1225 Toll Free Fax

104 John Kelly Drive,
Burnet, TX 78611
information@sw.stealthproducts.com
www.stealthproducts.com
facebook.com/StealthProducts
twitter.com/stealthproducts

Diego System Order Form



PO#			
Customer:			
Branch:		Date:	
Contact:			
E-mail:			
Phone:		Fax:	
Mark For:		Type:	<input type="checkbox"/> Quote <input type="checkbox"/> Order

Chair Info	Color: _____	Frame Width: _____	Frame Brand: _____	Frame Model: _____
Cane Size: <input type="checkbox"/> 3/4in <input type="checkbox"/> 7/8in <input type="checkbox"/> 1in				
Important Note: Some calculations on this order form are not automatically populated when an item is checked or selected, if a selection does not reflect a price change, move to next field by pressing tab which will allow the total cost to refresh. Prices and configurations subject to change without notice.				
Is this frame WC19 certified? <input type="checkbox"/> Yes <input type="checkbox"/> No				
[1] Warning: WC20 transit hardware requires: WC18 Restraint, WC19 Base and WC20 Seat/Back. Per RESNA WC-4:2012/ISO 16840. Click for more information.				

Step 1. Contoured Back [E2293]*	MSRP	WIDTH	HEIGHT	MOUNTING
<input type="checkbox"/> DGB Diego Contoured Back System with Coolcore® Cover and FlipLock™ Mounting Hardware				
See Back/Seat matrix table below Package Breakdown: DGMHFLB2 2 point FlipLock Universal Configuration Mounting Hardware DGBAL Back Cushion + Base + Cover				
Contoured Back System includes FlipLock hardware, contoured back cushion and base/wings. Foam package will affect final MSRP. The base MSRP includes Diego Back shell with cushion, Coolcore cover and package A foam, which consists of 1/2in HR3 contact foam and 1/4in base Closed Cell foam.				

Step 1.1 Back Foam Options	Step 2. Lateral Support [E0956/E1028]*	MSRP	MSRP
<input type="checkbox"/> Package A (included) 1/2in HR3 Contact Foam 1/4in Closed Cell Base Foam	Lateral Swing Away Hardware <input type="checkbox"/> DGLTSA-L Diego Lateral Swing Away Left \$280.00 <input type="checkbox"/> DGLTSA-R Diego Lateral Swing Away Right \$280.00	Lateral Fixed Hardware <input type="checkbox"/> DGLTFX-L Diego Lateral Fixed Left \$130.00 <input type="checkbox"/> DGLTFX-R Diego Lateral Fixed Right \$130.00	
<input type="checkbox"/> Package B (add \$95) 1/2in SunMate® Medium Contact Foam 1/2in HR3 Base Foam	Lateral Swing Away Pads (price shown per pair) <input type="checkbox"/> DGLPS33C Diego Lateral Pad for Swing Away 3x3in \$105.00 <input type="checkbox"/> DGLPS34C Diego Lateral Pad for Swing Away 3x4in \$105.00 <input type="checkbox"/> DGLPS44C Diego Lateral Pad for Swing Away 4x4in \$105.00	Lateral Fixed Pads (price shown per pair) <input type="checkbox"/> DGLPF33C Diego Lateral Pad for Fixed 3x3in \$105.00 <input type="checkbox"/> DGLPF34C Diego Lateral Pad for Fixed 3x4in \$105.00 <input type="checkbox"/> DGLPF44C Diego Lateral Pad for Fixed 4x4in \$105.00	

Step 3. Contoured Seat [E2294]*	MSRP	WIDTH	DEPTH	MOUNTING
<input type="checkbox"/> DGS Diego Contoured Seat System with Coolcore Cover				
See Back/Seat matrix table below. Package includes contoured cushion + seat shell + mounting hardware.				

Back & Seat Matrix	Step 4. Head Support [E0955/E1028]*	Step 5. i-Fit™ Pelvic & Posture Support																									
<table><tr><th>in (cm)</th><th>10 (25.40)</th><th>12 (30.48)</th><th>14 (35.56)</th><th>16 (40.64)</th></tr><tr><th>Width</th><th>Height / Depth</th><td></td><td></td><td></td></tr><tr><td>10 (25.40)</td><td>✓</td><td>✓</td><td>✗</td><td>✗</td></tr><tr><td>12 (30.48)</td><td>✗</td><td>✓</td><td>✓</td><td>✗</td></tr><tr><td>14 (35.56)</td><td>✗</td><td>✗</td><td>✓</td><td>✓</td></tr></table>	in (cm)	10 (25.40)	12 (30.48)	14 (35.56)	16 (40.64)	Width	Height / Depth				10 (25.40)	✓	✓	✗	✗	12 (30.48)	✗	✓	✓	✗	14 (35.56)	✗	✗	✓	✓	<input type="checkbox"/> CP06CCH51 \$459.00 Includes: CP06CCFPR Pediatric Comfort Plus 6in Pad \$151.94 HMO501 Multi-Axis Universal Pad Mount \$309.75 <input type="checkbox"/> CP08CCH51 \$459.00 Includes: CP08CCFPR Pediatric Comfort Plus 8in Pad \$151.94 HMO501 Multi-Axis Universal Pad Mount \$309.75	2 Point Pelvic Belts [E0978]* <input type="checkbox"/> IFPB5111 Push Button, 1in (2.54cm), Single Pull X-Small Pads \$86.00 <input type="checkbox"/> IFPB5112 Push Button, 1in (2.54cm), Single Pull Small Pads \$86.00 4 Point Pelvic Belts [E0978]* <input type="checkbox"/> IFPB5145 Push Button, 1in (2.54cm), Single Pull X-Small Pads \$130.00 <input type="checkbox"/> IFPB5146 Push Button, 1.5in (3.81cm), Single Pull Small Pads \$130.00 Posture Support [E0960]* <input type="checkbox"/> IFPS1101 Classic, Dynamic, 1in (2.54cm), Pull Down, X-Small \$142.00 <input type="checkbox"/> IFPS1102 Classic, Dynamic, 1in (2.54cm), Pull Down, Small \$142.00 <input type="checkbox"/> IFPS1201 Classic, Structure, 1in (2.54cm), Pull Down, X-Small \$142.00 <input type="checkbox"/> IFPS1202 Classic, Structure, 1in (2.54cm), Pull Down, Small \$142.00 <input type="checkbox"/> IFPS2101 Contoured, Dynamic, 1in (2.54cm), Pull Down, XS \$142.00 <input type="checkbox"/> IFPS2102 Contoured, Dynamic, 1in (2.54cm), Pull Down, S \$142.00 <input type="checkbox"/> IFPS2201 Contoured, Structure, 1in (2.54cm), Pull Down, XS \$142.00 <input type="checkbox"/> IFPS2202 Contoured, Structure, 1in (2.54cm), Pull Down, S \$142.00
in (cm)	10 (25.40)	12 (30.48)	14 (35.56)	16 (40.64)																							
Width	Height / Depth																										
10 (25.40)	✓	✓	✗	✗																							
12 (30.48)	✗	✓	✓	✗																							
14 (35.56)	✗	✗	✓	✓																							

i If you need a more customizable system not available with the Diego System consider our TFB (True Forming Back) Series, sizes start from 10in up to 20in wide/height. For more information call your Stealth Representative or one of our expert CSRs.

Order total based on products MSRP, any international or distributor discount is not applied, discounts will be applied to invoice. To inquire about any special discount you need to contact your CSR or local Stealth's representative.

ORDER TOTAL